

# **IMPACT**<sup>®</sup> **On-Board Aspirator SPECIFICATIONS**

**Controls:** OFF/On, Vacuum Adjust

**Indicators:** Vacuum Gauge, Dual Scale (0-760mmHg, 0-30" Hg)

**Collection Canister Capacity:** 1.1 Liter (1.15Qt.)

**Collection Canister Type:** Reusable Jar (may be used with or without disposable bag liners)

**Mounting Configuration:** Wall or Rail

**Warranty:** One (1) Year, Limited Parts and Labor

**Included Accessories and Installation Component Parts:**

**Tubing:** 1 ea. 5' x 3/8" ID & 10" x 3/8" ID

**Disposable Bag Liners:** 1 dz.

**Fittings:** 1 ea. "Y"

2 ea. Suction Tips

1 ea. Quick-Disconnect (319K only)

1 ea. Hose Barb, 3/8" ID Tube (329K only)

1 ea. Airway; 60 mm, 80 mm, 100 mm, (329K only)

1 ea. Catheters; 14 French, 18 French (329K only)

1 ea. Rinse Bottle (329K only)

**Manual:** 1 ea. Model 319/319k

**Warranty Card:** 1 ea.

**Miscellaneous:** 5/pkg Hydrophobic Membrane Filter (optional)

**Note:** Impact Instrumentation Inc., reserves the right to change these specifications without prior notification.

p/n 929-0319-00 Rev A      Updated 10/11

**Impact Instrumentation, Inc. - P.O. Box 508 - West Caldwell, NJ 07007-0508 Tel: 973.882.1212 Fax: 973.882.4993**