

IMPACT[®] 324/ 324K On-Board Aspirator SPECIFICATIONS

Vacuum: Adjustable, 0 to 550 mmHg (0 to 22" Hg)

Free Airflow: 30 liters per minute at patient end of suction tube (min)

Controls: Power OFF/On, Vacuum Adjust

Indicators: Vacuum Gauge, Dual Scale, 0-760mmHg (0-30" Hg) Lamp, Illuminates when powered "ON"

Collection Canister Type/Capacity: Disposable, 1200ml (includes bacterial filter)

Mounting Configuration: Action-area wall, collection-canister-to-the-left or collection-canister-to-the-right (as specified)

EMI/RFI Certification: To Mil-Std-461/462 (Model 324K only)

Warranty: Two (2) Years, Limited Parts and Labor

Included Accessories and Installation Component Parts:

- 1 ea. 820-0018-00 Suction hose, Clear
- 1 ea. 540-0163-00 Exhaust Hose, 5' x 3/8" ID, Black
- 2 ea. 410-0004-00 Collection Canister/Lid, Disposable, 1200ml
- 1 ea. 906-0324-02 Manual, Operation, Model 324/324K
- 1 ea. 940-0001-00 Warranty Registration Card
- 1 ea. 820-0031-00 Yankauer Suction Tip with Thump Port (Model 324K only)
- 1 ea. 820-0024-00 Rinse Bottle, 8 oz. (Model 324K only)
- 1 ea. 703-0324-09 Assembly, Mounting Bracket (includes hardware items below)
- 4 ea. 358-1032-08 Screw, Phillips, Pan Head, 10-32 x 1/2
- 4 ea. 376-0020-00 Washer, Lock Internal Tooth, #10

Note: Impact Instrumentation Inc., reserves the right to change these specifications without prior notification.

p/n 929-0324-00 Rev A Updated 10/11