

# IMPACT 305 VAC-PAK SPECIFICATIONS:

VACUUM RANGE:	0-550 mmHg (0-22" Hg)
FREE AIR FLOW:	Minimum 31 Liters Per Minute (LPM)
TEMPERATURE OPERATING RANGE:	-60°C to 60°C (-76°F to 140°F)
INPUTS:	External AC/DC Charger Adapter External 12 VDC
CONTROLS:	Power Switch- OFF/CHG-EXT-INT
DISPLAYS:	Charge Indicator Lamp
COLLECTION CAPACITY:	500cc (.50 Liter)

## POWER

EXTERNAL 117 VAC:	Recharge only with AC/DC Charger Adapter: PICB 133, 810-0001-00, Shock Protection: Type B, Class 2
EXTERNAL 117 VAC:	Continuous operation or recharge (with optional AC charge rectifier)
EXTERNAL 12 VDC:	Continuous operation or recharge
INTERNAL OPERATING TIME:	40 minutes
INTERNAL RECHARGE TIME:	14-16 hours

## ACCESSORIES

STANDARD:	DC CHARGER CABLE, AC/DC CHARGER ADAPTER, 6 FT. CLEAR SUCTION TUBING, RINSE BOTTLE, 2 SUCTION TIPS, "Y" FITTING, 2 CATHETERS, 9" CLEAR HOSE
OPTIONAL:	CATHETERS (SPECIFY SIZE & STYLE), OROPHARYNGEAL AIRWAYS, "T" FITTINGS, AC CHARGER RECTIFIER (SPECIFY 117 VAC, 60 HZ - OR - 220-240 VAC, 50 HZ)
CASE MATERIAL:	Nylon, Cordura, Water-Resistant
CASE SIZE:	16.5 cm W x 16.5 cm D x 28 cm H (6.5" W x 6.5" D x 11" H)
WEIGHT:	4.1 kg 9 lbs
WARRANTY:	Limited, 1-year (see LIMITED WARRANTY statement)

Note: Impact Instrumentation Inc., reserves the right to change these specifications without prior notification.

Rev A (08/06)